

PHCOM030

PHCOM030 Caretex® Commode Pan Liner 1700ml

caretex

PRODUCT INFORMATION

ltem Number	РНСОМозо
Product Colour	Grey
Volume	1700ml

FEATURES

Designed to be used with commode chairs

Easy to store as each product is stackable and separate

Disposable product for increased hygiene and reduced risk of cross contamination

Made of recyclable materials

Rack available for use with this product

Steady when standing due to flat bottom surface.

Strong sides able to support weight of contents when full

SECTORS

Medical, Vet and Dental • Pharmaceutical & Laboratory

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. You are advised to retain the packaging for reference. Always inspect product prior to use and do not use if damaged.

PACKING

200 pieces per case 10 cases per pallet

DISPOSAL INSTRUCTIONS

Dispose of responsibly and in line with local guidance.





PHCOM030

PHCOM030 Caretex® Commode Pan Liner 1700ml



Size	Outer Packaging Barcode
One Size	5024951012825

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

Polyco Healthline, South Fen Business Park, South Fen Road, Bourne, Lincolnshire, PE10 ODN Tel: 0333 320 8550 • E: sales@polycohealthline.com • www.polycohealthline.com